


ADD6 JANUARY 2026

AMORPHOUS PHARMACEUTICALS

CHRIS BENMORE

Senior Physicist
Advanced Photon Source, Argonne National Laboratory, IL
Arizona State University, Tempe, AZ
University of Chicago, Chicago, IL




U.S. DEPARTMENT
of **ENERGY**

Argonne National Laboratory is a
U.S. Department of Energy laboratory
managed by UChicago Argonne, LLC.

Argonne 
NATIONAL LABORATORY

Advanced
Photon Source

ILL/ESRF PARTIAL STRUCTURE FACTORS WORKSHOP 1992



COLLABORATORS & ACKNOWLEDGEMENTS



Pam Smith & Steve Byrn



Jeffery Yarger




Rick Weber & Stephen Wilke



Funding :

- National Institutes of Health, R44GM117701.
- U.S. DOE Office of Science, DE-AC02-06CH11357
- Argonne, LDRD



**SPS-XRPD
WORKSHOP**

2018 Purdue University, USA
2019 PSI, Switzerland
2023 Purdue University, USA
2024 Novartis, Basel, Switzerland
2025 APS, Argonne, USA
2026 Novartis, Basel, Switzerland



EXCELSUS
STRUCTURAL SOLUTIONS




Fabia Gozzo, Thibaud Stoll & Mathilde Reinle-Schmitt




Arnaud Grandeury

SPS-XRPD 2025 attendees



SOLUBILITY, BIOAVAILABILITY, & FORMULATION





U.S. DEPARTMENT
of ENERGY

Argonne National Laboratory is a
U.S. Department of Energy laboratory
managed by UChicago Argonne, LLC.

Argonne 
NATIONAL LABORATORY

Advanced
Photon Source


RITONAVIR



Law *et al.*, J. Pharm. Sci. 93 (2004) 563

- Originally dispensed in 1996 in crystalline form I (monoclinic)
- HIV deaths in US fell by 64% following FDA approval of Ritonavir
- Two years later a lower energy more stable form II appeared (orthorhombic)
- Form II was much less soluble and had lower bioavailability
- Caused the removal of formulation from the market (cost Abbvie ~US\$240M)
- Re-formulated as amorphous form in 2000 with higher solubility.

BIOPHARMACEUTICS CLASSIFICATION SYSTEM





- Low aqueous solubility is the major problem encountered with formulation development.
- More than 40% of new drugs developed in the pharmaceutical industry are practically insoluble in water.
- Solubility, permeability, polymorphism and stability are key issues in pharmaceutical formulations
- Class II pharmaceuticals: Posaconazole, Itraconazole, Carbamazepine, Nifedipine, Felodipine, Ketoprofen

ACOUSTIC LEVITATION



A LABORATORY MODEL FOR SPRAY DRYING

- Spray drying produces dry powder from liquid by rapid drying with hot gas
- Spray nozzle disperses droplets in typically in 100-300 micron size range
- Commonly used method for producing amorphous pharmaceuticals.
- Acoustic Levitator operates at 22 kHz
- Nodes are 7mm apart, so sample size is limited up to 3mm in diameter
- Surface tension ultimately determines droplet diameter



C.J. Benmore and J.K.R. Weber *Containerless synthesis of amorphous and nanophase organic materials*. U.S. Patent 9327264B2 (2016)



AMORPHOUS SOLID DISPERSIONS PREPARED BY DIFFERENT METHODS



- Amorphous Solid Dispersions are used to increase the solubility of oral medicines by stabilizing the more soluble amorphous phase
- Ketoprofen (KTP) – Polyvinylpyrrolidone (PVP) mixtures made by different methods – how can PDF help ?


LevPM - premix levitation
LevCS - co-spray levitation
MQ - melt quenching
RV - rotary evaporation
SD - spray drying

PHASE SEPARATION IN KTP-PVP ASD'S



- Variability in the extent of molecular mixing using different amorphization routes - a direct indication of their relative stability against crystallization.
- MQ, SD & LevPM solutions exhibit some degree of molecular mixing. LevCS close to KTP/PVP physical mixtures

EPSR MODELING OF ACTIVE PHARMACEUTICAL INGREDIENTS (API)



U.S. DEPARTMENT
of ENERGY

Argonne National Laboratory is a
U.S. Department of Energy laboratory
managed by UChicago Argonne, LLC.


Argonne
NATIONAL LABORATORY



Advanced
Photon Source

HIGH ENERGY X-RAY DIFFRACTION

- 6-ID-D at the Advanced Photon Source has a High Energy Bent Double Laue Monochromator
- Si 311 (50-95keV) Si 331 (70-125keV)
- Undulator 13th harmonic is at 90.88keV @ 9mm gap
- Calibrate to Bi K-edge 90.54 keV
- Data reduction for liquids problem: not all x-ray PDF software are the same.




Quantum Beam Sci. 2023, 7, 20.

Review

Review of Current Software for Analyzing Total X-ray Scattering Data from Liquids

Leighanne C. Gallington ¹, Stephen K. Wilke ¹, Shinji Kohara ² and Chris J. Benmore ^{1,*}



EMPIRICAL POTENTIAL STRUCTURE REFINEMENT

- EPSR is a RMC approach: a fitting tool analogous to Rietveld refinement for liquids
 - Uses a reference LJ potential & atomic charges PLUS an Empirical Potential
 - Empirical Potential altered until model and diffraction (PDF) data agree


$$U_{\alpha\beta}(r) = 4\epsilon_{\alpha\beta} \left[\left(\frac{\sigma_{\alpha\beta}}{r} \right)^{12} - \left(\frac{\sigma_{\alpha\beta}}{r} \right)^6 \right] + \frac{q_{\alpha}q_{\beta}}{4\pi\epsilon_0 r} + c_{\alpha\beta} \exp \left[\frac{1}{\gamma} (r_{\alpha\beta} - r) \right]$$

Intermolecular
potential energy

LJ potential


Coulomb

Empirical potential



MODELING PROCEDURE - EPSR

- Use LJ starting parameters from Optimized Potentials for Liquid Simulations (OPLS). Listed by chemical group (W.L. Jorgensen)
- Partial Charges taken from MD simulations or DFT calculations
- Define torsional angle rotations of specific bonds from NMR or polymorphs of specific chemical groups
- General strategy – Test a range of models. Focus on simplest models that provide a “good” fit to $S(Q)$.
- Try different starting molecular conformations
- Typically 60-100 molecules in a cubic $\sim 30\text{\AA}^3$ box
- Collect over 10,000 configurations – typically 12 hours



RIGID, ROD-LIKE MOLECULES

Define **Posaconazole** molecule in EPSR




Atom types (chem groups): 26

Angles (3 atoms): 90

Dihedral angles (4 atoms): 249

Rotations specified: 13 tried – rigid works well



Standard LJ parameters: changing has little effect





A SAXS-WAXS study of the
endothermic transitions in amorphous
and supercooled liquid **Itraconazole**.

C.J. Benmore et al.

Thermochimica Acta 644 (2016) 1



RIGID MOLECULES - CARBAMAZEPINE




- Carbamazepine (CBZ) used to treat epilepsy

- Liquid & Amorphous CBZ has a much higher oral solubility than crystal forms

- Stability is an issue. In 1998, 70 million tablets of CBZ were withdrawn – dihydrate formed


- N1-O1 Hydrogen bonding. C3-C3 ring interactions

- Liquid and glass models show structures not present in the crystalline forms that vary with temperature





CARBAMZEPINE



CBZ HYDROGEN BONDING & AROMATIC INTERACTIONS






Double HB dimers in crystal → mainly single HB bond in liquid/glass



More aromatic ring interactions in liquid than glass




More HB trimers and tetramers in liquid than glass





ROTATIONS MATTER – CBD (CANNABIS)

NMR shows rigid orientation between two rings 6-4
This is an essential constraint to fit diffraction data





- Systematic inclusion of rotations to observe effect on EPSR fit
- 8-15 has two different conformations in crystal
- Need to allow 12-19 rotation to reproduce FSDP
- Sparse HB from H10. H13 partially blocked by side group




NIFEDIPINE & FELOPIDINE

- Nifedipine (NIF) and Felodipine (FEL) are both calcium channel blockers - strong cardiac medicines for hypertension
- NIF has a nitro group. FEL contains a dichlorophenyl group.
- NIF crystallizes faster than FEL. Questions regarding relative strengths of hydrogen bonds.
- Partial charges were assigned based on MD simulations. Gupta & Kotheke, *J. Biosci.* 22 (1997) 177.
- Charge balance on each chemical group to probe chemical group swap.
- Full rotations and partial rotations (wobble) introduced (no dihedrals defined but no rotations enabled)




NIF POLYMORPHS


GLASSY NIFEDIPINE





- CCD show NIF has six polymorphs. 10 conformers
- α -phase is the most stable and grows $>120^\circ\text{C}$
- $<42^\circ\text{C}$ liquid NIF turns to glass and its crystallizes to the β phase
- Models with α & β form conformations



α & β NIF MODELS




- Better fits $<6\text{\AA}$ for β -conformation -consistent with DFT - predicts most crystal conformers have similar energies.
- Weak inter-molecular hydrogen bonds of the type $\text{N-H}\cdots\text{O}$ in both the α and β NIF polymorphs.




FELODIPINE (FEL) POLYMORPHS



- FEL is a chiral drug with two enantiomers: R- & S-
- The R-enantiomer of FEL is more likely to cause adverse reactions than the S-enantiomer
- Form I is the most stable. Form II is metastable, has the highest solubility & dissolution rate. Both are racemic



R-enantiomer (form I)




S-enantiomer (form II)





- Five polymorphs
- Use racemic mixture of R & S based on form I conformation
- $T_m \sim 417K$ & $T_g = 318K$


RACEMIC FEL MODEL



- Disorder in side chains between form I and form II reflected between enantiomers in 2nd peak in $D(r)$
- Previous X-ray PDF and DFT suggest extra intra-molecular alkoxy–Cl correlations at 3.0 \AA . Geddes *et al. Chem. Commun.*, 2019,55, 13346.
- Our EPSR shows 3.0 \AA distance corresponds primarily to *intra-molecular Cl-Cl pairs*.
- Slight shift corresponds to 4° change in $\angle\text{CCCl}$ angle within the dichlorophenyl ring




BONDING IN GLASSY NIF & FEL




- Liquid and glassy NIF & FEL both have N1H-O2 hydrogen bonds
- **NIF also has similar number of N1H-O3 (nitro) hydrogen bonds – an interaction not observed in the crystal forms**
- FT-IR has previously suggested hydrogen bonding in the amorphous form of FEL is stronger than in crystal. Marsac *et al. Pharm Res.* 23 (2006) 2306.
- Not supported by EPSR
- Based on our EPSR models we suggest a peak re-assignment to “free” (non-hydrogen bonded) N-H stretching groups

AMORPHOUS SOLID DISPERSIONS






POLYMER “EXCIPIENT”

- PVP is a **hydrophilic** polymer commonly used as an excipient (inactive ingredient) in formulations
- In Amorphous Solid Dispersions, excipients inhibit crystallization (increase shelf life) and promote stability. Effectiveness depends on T_g , molecular mobility & API-polymer bonds
- Hydrogen bonding interactions between the API & excipient are often the main cause of the dissolution enhancement in the bloodstream



WATER IN PVP







- Water uptake depends on temperature and humidity. Don't store in the bathroom!
- EPSR shows linear relations between the water-water and water-polymer CN's versus water content in PVP.
- A stronger preference for water-water hydrogen bonding over water-polymer hydrogen bonding is found.

EPSR SENSITIVITY TO INPUT PARAMETERS

KETOPROFEN-PVP


- Assess EPSR sensitivity to:
 - Intramolecular bond rotations
 - PVP molecule length
 - PVP tacticity

Simulation	Intramolecular rotations [*]	PVP length (# monomers)	PVP tacticity
A	none	5	syndio
B	2, 5, 8, 11, 12	5	syndio
C	1-5, 6, 8, 10-13	5	syndio
D	none	3	syndio
E	none	10	syndio
F	none	5	iso





- Six simulations used to probe EPSR's sensitivity to different molecular inputs
- Includes 9.4 wt.% water

GOODNESS-OF-FIT COMPARISONS



- R -factors indicate better goodness-of-fit in real space
- A, B & C explore **bond rotations**:
 - A has rigid molecules – good fit
 - B & C have increasing numbers of bonds rotations – worse at high- Q (intra-molecular problems in PDF)
- A(5), D(3) & E(10) explored the effect of **PVP molecule length** – all rigid.
 - R -factors decrease as the polymer length decreases.
- Pre-peak $Q=0.85 \text{\AA}^{-1}$ from the PVP chain.
- First diffraction peak in $S(Q)$ is best matched to (D) 3-monomer PVP
- F used rigid, 5-monomer, **isotactic PVP**, which yielded the worst match to X-ray data.



BOND DIHEDRAL ANGLE DISTRIBUTIONS



- EPSR allows a few degrees of bond rotation to simulate thermal motion.
- narrow dihedral angle distributions in A.
- Bond rotations enabled in B & C but interatomic potentials prevent intra-molecular overlap.
- When bonds rotate, either a few distinct angles are adopted due to steric hinderance e.g. purple 5 in B, or a nearly random distribution occurs e.g. blue 10 in C.

HYDROGEN BONDING INVOLVING KTP

- Hydrogen bonds are distinguished by the second molecule : KTP molecule, PVP, or water.
- (A) KTP molecule has 1.78 hydrogen bonds –majority are water bonds.
- All six simulations similar except (D) 3-PVP monomers



- The number of KTP-KTP hydrogen bonds increases (0.20 to 0.69) with drug loading.
- KTP-PVP hydrogen bonding is similar with composition – indicates reported instability of high drug loading ASDs likely due to a decrease in PVP's steric hinderance of the KTP molecules

FUTURE DEVELOPMENTS

SRN

Synchrotron Radiation News
July/August 2022 • Vol. 35, No. 4

 Taylor & Francis
Taylor & Francis Group



Machine Learning



U.S. DEPARTMENT
of **ENERGY**

Argonne National Laboratory is a
U.S. Department of Energy laboratory
managed by UChicago Argonne, LLC.

Argonne 
NATIONAL LABORATORY

**Advanced
Photon Source**


PROS & CONS

- EPSR is a powerful technique for investigating disordered structures – input knowledge from known polymorphs, NMR, vibrational spectroscopy and DFT calculations
- EPSR is not unique but ideal for testing different molecular conformation/bonding scenarios
- Main drawback is it's classical by nature- no quantum mechanics are included
- Both scientific and legal implications to this work. How to patent an Amorphous Solid Dispersion ?




- Ibuprofen is a flexible molecule with two polymorphs. Form I was discovered in 1974 and form II in 2010
- “Blind test” predicted form II crystal structure (6th)
- Francia *et al.* (2021) studied 555 possible crystal structures, optimizing both the molecular conformation and crystal structure, and found 205 viable structures!

FUTURE DIRECTIONS – MACHINE LEARNING



- *Deciphering diffuse scattering with machine learning...* G. Sivaraman & C.J. Benmore *J. Phys. Condens. Matt.* 36 (2024) 381501

- Experimentally Driven Automated Machine-Learned Interatomic Potential for a Refractory Oxide. Sivaraman... C.J. Benmore. *Phys. Rev. Lett.* 126 (2001) 156002



- Tested different combinations of torsional angles of rotation, τ_n , within the Ibuprofen molecule.
- Quantum chemical calculations based on machine learning interatomic potentials (Rowan software) were performed on the Ibuprofen molecule and found six likely conformers – EPSR input
- Rowan tools do not possess long range forces needed for polymer reactivity, complex melts.






Argonne
NATIONAL LABORATORY




Advanced
Photon Source


High Energy Bent Double Laue Monochromator





2mm CVD
diamond
4mm Al filters



1st Crystal in
liquid bath of
Ga-In-Sn





2nd Crystal




6-ID-B,C
Bragg Monochromator

6-ID-D
Laue-Laue Monochromator



NIFEDIPINE (NIF) POLYMORPHS


- CCD shows NIF has six polymorphs. α -phase is the most stable and grows $>120^\circ\text{C}$
- $<42^\circ\text{C}$ liquid NIF turns to glass and its crystallizes to the β phase




FELODIPINE (FEL) POLYMORPHS

- FEL is a chiral drug with two enantiomers: R- & S-
- The R-enantiomer of FEL is more likely to cause adverse reactions than the S-enantiomer
- Interchange of the methyl and ethyl chains inverts the chirality of the molecule
- Form I is the most stable and racemic
- Form II is metastable, has the highest solubility & dissolution rate


R-enantiomer




S-enantiomer




Form I
DONTIJ




Form II
DONTIJ01



Form III
DONTIJ02



Form IV
DONTIJ11



Form V
DONTIJ06

